

[www.zhichenginstrument.com](http://www.zhichenginstrument.com)

**Model: ZHWY-3312** *New for 2012!*

**Extra Large Capacity Orbital Shaker**

(Note: The actual product looks much better than the below photo.)



**Maximum  
Capacity**

**Three Layers  
Heavy Duty  
Support 5000mlx36  
flasks**

## Technical Specifications

<b>Model</b>	<b>ZHWY-3312</b>
<b>Item Name</b>	<b>Orbital Shaker</b>
Controller	P.I.D microprocessor
Shaking Speed	30-180 rpm ( $\pm 1$ rpm)
Shaking Mode	Orbital
Diameter of Orbit	50mm
Timer	0-500 hours
Setting	Digital
Display	LCD
Capacity	250mlx288 or 500mlx180 or 1000mlx120 or 2000mlx72 3000mlx48 or 5000mlx36
Platform Size	780x620 mm
Platform included	6
Exterior Dimension	1540x980x1293 mm
Voltage	220V 50/60Hz
Power	1000W
Weight	630Kg
Warranty	1 year