

www.zhichenginstrument.com

Model: ZHWY-A2005 Shaking Incubator


New!

- *Space Saving*
- *Stackable*
- *Double Chambers of independent control*
- *Slide-out platforms*



Note: The actual product looks much better than the above photo.

Control Panel



Technical Data & Specifications

Model	ZHWY-2005A incubator shaker
Shaking Mode	Orbital
Refrigeration	CFC-free, R134
Temperature Range	Ambient-20°C to 60°C (or 4-60°C at ambient 25°C)
Temperature Set Division	±0.1°C
Temperature Uniformity	±1°C
Heat Distribution	Forced-air circulation
Controller	P.I.D Microprocessor
Setting	Digital
Display	LCD
Diameter of Orbit	50mm
Shaking Speed	30-300 rpm (±1 rpm)
Timer	1 minute to 500 hours
Capacity/each chamber	100mlx50 or 250mlx32, 500mlx26 or 1000mlx14 2000mlx8
Platform Size	800x430 mm
Platform included	2
Chamber Size/each	920x532x373mm (182L)/each
Motor Power	2400W
Electrical Requirement	AC 220V 50/60Hz
Weight	440Kg
External Dimension	110x91x133cm
Warranty	1 year